

UNIVERSITY OF
BIRMINGHAM

U


The Wear of Ceramic Total Disc Replacements

B

Naomi Green,
Athena Swan Event
12th November 2014

Background to Research


- Back pain is a debilitating condition
- Costs UK economy billions of pounds a year in treatment and lost working days



Anatomy of a Spinal Segment [1]



Healthy Disc [1]



Extreme Degeneration [1]

[1]: S Kurtz, A Edidin, Spine Technology Handbook

UNIVERSITY OF
BIRMINGHAM

Total Disc Replacements

- Total Disc Replacement (TDR) possible treatment
- Intervertebral disc is removed and replaced with TDR
- Issues with devices and debate over clinical outcomes
- Ceramic materials may improve devices



Synthes Spine ProDisc II: Total Disc Replacement [1]



Medtronic MAVERICK: Metal on Metal TDR [1]



Scient'x Discocerv: Ceramic on Ceramic TDR [2]

[1]: S Kurtz, A Edidin, Spine Technology Handbook,

[2] http://www.scientx.com/product_nonfusion_discocerv.php 2006

UNIVERSITY OF
BIRMINGHAM


Aims

- Overall aims:
 - Investigate the use of ceramics in cervical Total Disc Replacements and their wear behaviour
 - Manufacture ceramic implants using novel grinding methods
 - Develop a new method for measuring low wear rates for implants

UNIVERSITY OF
BIRMINGHAM

Implant Manufacture

- Universitas 21 Scholarship to study at University College Dublin, Ireland
- Collaborative study of the grinding of ceramic TDRs using novel diamond tool to achieve desired surface finish



Manufacture of Implant




Upper and Lower Alumina Implants in Fixtures

UNIVERSITY OF
BIRMINGHAM

Wear and Friction Testing

- Major mode of failure for an implant is wear
- TDRs are tested on spinal simulators to the international standard ISO 18192-1 [3] to determine wear rates



Bose SDWS-1 Spine Simulator

[3] ISO 18192-1:2011 Implants for surgery — Wear of total intervertebral spinal disc prostheses Part 1

UNIVERSITY OF
BIRMINGHAM


Measuring Wear Rates

- Normal method for measuring wear rates of implants is the gravimetric method
- Change in mass is used to determine volume loss
- When wear rates are low it is difficult to measure the change in mass of the sample accurately
- Alumina cervical implants are expected to have low wear rates
- Alternative method is required
- Use non-contacting optical method, Alicona Infinite Focus, to scan surface of the implant
- Can measure form and surface roughness parameters in one scan

UNIVERSITY OF
BIRMINGHAM

Implant Wear Measurement

- Form and roughness measurements taken before start of wear test and at every 0.5 million cycles
- Matlab code has been developed which finds volume difference between the unworn and worn surfaces




Images of the Functional Surfaces Before and After Testing

UNIVERSITY OF
BIRMINGHAM

Friction Measurements

- Short 100 cycle tests
- Measure frictional torque against frequency
- Produce Stribeck Curves
- Gradient of the line gives information about lubrication regime



Example Stribeck Curve

UNIVERSITY OF
BIRMINGHAM

Future Work

- Complete final scans of implant
- Analyse data
 - Calculate wear rates
 - Compare gravimetric method with optical method
 - Surface roughness changes with wear
- Publish journal papers
- Write thesis
- Graduate !?!

UNIVERSITY OF
BIRMINGHAM